ARTIFICIAL INSEMINATION - INSEMINATION

MULTI-DOSE KIT

Backpack for the modulated multi dose of semen extender. Capability to contain bags of 3 or 5 lt of liquid. Ideal for multiple inseminations using the dosing syringe. Available big syringe with pipes and fittings.

CODE	DESCRIPTION
1907002	Backpack for multi dose
3944080	Big syringe with pipes and fittings
1907020B	Disposable bag 3 lt (pack 20 pcs.)
1907113B	Plug for Peristaltic Pump (pack 20 pcs.)
1907115	Plastic clip for stops liquid flow




SPRING FOR AUTO-INSEMINATION

Spring for hanging aluminium bag or insemination tube.



PLASTIC HOLDER

Plastic holder with iron bar for attachment of the probe. As an alternative to the use of the bag, allows to stimulate the sow during artificial insemination, even in the absence of the operator.



BOARMATE SPRAY

Spray with boar odour to detects sows heat, or use as an extra stimulant during insemination. Spray maximum for 2 seconds near the nose of the sow.



ARTIFICIAL INSEMINATION - INSEMINATION


CATHETERS FOR A.I.

COD. 1906009

COD. 1906010

Catheters for artificial insemination, with reclosable plug. Available in different packagings and colors.

CODE	DESCRIPTION
1906007	Yellow sponge catheter with plug, with connector 2, neutral box (pack 5 pcs.)
1906008	Yellow sponge catheter with plug, easy opening, neutral box (single pack)
1906002	White sponge spiral catheter (pack 5 pcs.)
1906009	Green sponge catheter for young sows, with connector 2, (pack 5 pcs.)
1906010	Green sponge catheter for young sows, with connector 2, (single pack)
1906014	Blue flexible connecting tube (pack 10 pcs.)
1906018	White catheter in anti-reflux silicone (single pack)




INTRAUTERINE CATHETERS AND TUBES

Catheters and tubes for the intrauterine artificial insemination, without reclosable plug. Available in two different models.

CODE	DESCRIPTION
1906006	Green intrauterine sponge catheter (single pack)
1906013	Green intrauterine tube, with two holes in the side (single pack)
1906020	Green intrauterine sponge catheter, with three holes in the side and cap (single pack)
1906021	Green intrauterine sponge catheter, anti-reflux system with cap (single pack)



ARTIFICIAL INSEMINATION - INSEMINATION



ARTIFICIAL INSEMINATION -INSEMINATION

CLEANING WIPES

Wipes to clean the udders and the vulva of the sow before the artificial insemination. Useful for the stimulation of milk production and to reduce the time of milking. Prevents ascending contaminations in the reproductive apparatus. Suitable also for the cleaning of the operator's hands. Only for external use. Bucket of 1000 wipes. Measures: 20 x 20 cm.

CODE	DESCRIPTION
4980030	Cleaning wipes
4980031	Wipes refill



COD. 4980030

VULVA CLEANSER

No spermicidal cleansing solution that allows for a fast cleaning of the vulva before the insemination. Apply a few sprays on the part to be cleaned and dry up with a clean paper cloth. In most cases this operation avoids the occurence of ascending infections in the sow reproductive apparatus. Available in two different packages.

CODE	DESCRIPTION
1907078	Vulva cleanser in bottle 1lt
1907079	Vulva cleanser in tank 5 It



LUBRICANTS

No spermatoxic gel to lubricate the catheter before intserting it. It also has a bacteriastatic power, therefore it reduces the bacterial charge in the reproductive apparatus. It can be used by inserting 10 - 15 ml in the cervix after the insemi-nation, thus reducing the risk of seminal material rejection. Available also lubricant oil, cheap and easy to use. Available in different packages.

CODE	DESCRIPTION
3946063	Lubricant gel in bottle 250 ml
3946062	Lubricant gel in Bag in box 5 lt
1907077	Vaseline oil in bottle 1 kg
1907076B	Vaseline oil in tank 5 kg



COD. 3946063

COD. 3946062